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
Place a clean microscope slide in a Petri dish. Place 5 drops of a 0.5 M glucose solution on the microscope slide. Add 25 drops of an active silver ion solution to the glucose solution. Wait 3 minutes while the solution darkens and a grayish precipitate forms. A silver mirror is also forming on the slide, though it may be obscured by the precipitate. Use water from a wash bottle to wash off the precipitate and reveal the silver mirror

Using a gloved hand to avoid contact with the silver ion solution, remove the slide from the Petri dish, rinse the silver mirror with water and let dry. (For faster drying a second rinse with acetone may be used.)


Add some water dropwise to the slide. How attracted are the water drops to the silver?

In a fume hood, cover the silver with a few drops of a long chain alkanethiol solution in ethanol.

Allow the ethanol to evaporate in a fume hood, leaving behind an alkanethiol monolayer with the sulfur atoms bound to the silver and the hydrocarbon tails pointing away. This effectively coats the surface with hydrocarbons. Add water dropwise as before. How attracted are the water drops to the surface? (Like attracts like.) Do water drops on the monolayer coated surface spread out or bead up?

Is the contact angle greater or less than before the alkanethiol was added? Is the water attracted more to the plain glass, to the silver, or to the alkanethiol monolayer-coated silver?

Materials

0.5 M glucose	Dissolve 0.90 g glucose in 10 mL of water.
0.8 M KOH	Dissolve 0.45 g KOH in 10 mL of water.
0.1 M silver nitrate	Dissolve 0.17 g AgNO ₃ in 10 mL of water.
15 M concentrated ammonia	
active silver ion solution	Add concentrated ammonia dropwise to 10 mL of 0.1 M silver nitrate solution until the initial precipitate just dissolves. Add 5 mL of 0.8 M KOH solution; a dark precipitate will form. Add more ammonia dropwise until the precipitate just redissolves. This "active silver" solution should be used within an hour of preparation. To avoid the formation of explosive silver nitride, discard any remaining active solution by washing down the drain with plenty of water.
Alkanethiol solution	Add a few drops of a long-chain alkanethiol, such as decanethiol or hexadecanethiol, to 20 mL of ethanol. Keep this solution in a fume hood because of the stench. When finished, dispose of this solution by adding about 5 mL of household bleach. Let stand for several minutes then wash the solution down the sink.

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
8 Ferrofluids

What are ferrofluids ?¹

Liquid magnets do not exist, since all magnetic substances lose their ferromagnetic properties above a certain temperature, the Curie temperature, which is always lower than the melting temperature. However, magnetic liquids do exist in the form of so-called ferrofluids; they are black liquids that can be manipulated by an external magnetic field gradient which attracts the liquid (fig 1). Ferrofluids consist of small magnetic colloids (which are particles of the order of 1 nm – 1 μm) in a non-magnetic solvent [1].

Colloidal particles are much larger than molecules and sediment under the influence of gravity. However, they are small enough to perform Brownian motion (due to the random thermal motion of the solvent molecules of the suspending fluid), which counteracts the sedimentation. Colloids usually attract each other due to Van der Waals forces, that have to be compensated by repulsive forces (e.g. charge or steric repulsion by surfactant layers) to obtain a stable colloidal dispersion, in which particles freely translate and rotate. In ferrofluids, as a special kind of colloidal dispersions, additional dipolar interaction forces between particles are present, that are direction dependent but attractive on average [1].

Figure 1. Example of an iron oxide dispersion in zero-field (left) and exhibiting the equilibrium shape in a field gradient typical for ferrofluids (right).

¹ Karen Butter, Dissertation Utrecht 2003, „Iron(oxide) ferrofluids: synthesis, structure and catalysis“; <http://www.library.uu.nl/digiarchief/dip/diss/2003-1027-144230/full.pdf>

Applications of ferrofluids²

Applications of magnetic fluids are numerous, and started around 1970. Many applications are based on the special property of ferrofluids, that they are attracted by magnetic field gradients while preserving their liquid character. They are used e.g. as liquid seals in rotary shafts, but also in computer hard disk drives, for prevention against contamination [2, 3, 5].

The same property is used in loudspeaker coils [5], that consist of a cylindrical voice coil fitted to a cylindrical permanent magnet with a small gap that allows the coil to move. The heat developed in the voice coil can be dissipated more readily when the air in the gap is replaced by a liquid. An ordinary liquid would drip out, but a ferrofluid is retained by the magnetic field already present.

Other applications of ferrofluids are the use as contrasting agent for MRI scans or as magnetic inks for jet printing [2], which are subjected to magnetic field gradients to deflect and position the ink droplets on the paper. Here, the magnetic properties of the dried ink are not relevant, but there are special constraints such as a quick drying time on the paper.

For magnetic reading of letters and numbers, e.g. on cheques, the magnetic properties of the dried ink are of paramount importance; the paper is first subjected to a magnetic field, before it is fed into the magnetic reader, so the dried ink must have a suitably high remanence for the letters to be read. Nowadays, research is performed to influence the color of magnetic inks by building fluorescent dyes into the particles [6].

Magnetic fluids are also used in medical applications, for instance in magnetically controlled drug targeting, based on the binding of e.g. established anticancer drugs with ferrofluids that concentrate the drug in the area of interest (a tumor) by means of magnetic fields [4].

Other potential applications of ferrofluids are their use as optical devices (e.g. shutters), exploiting the change in transmittance of light upon alignment of the particles in a magnetic field. Here, small magnetic particles are needed to decrease the response time of the shutter [3].

A relatively new potential application is the use of a lattice of self-assembled monodisperse magnetic colloids in computers for data storage [7,8]. However, new methods for writing and reading of data have to be developed first, to be able to use this in practice.

Another very interesting application would be the use of magnetic colloids attached to catalysts, which could be separated after their use by means of magnetic field gradients. This application would then combine the advantages of homogeneous catalysis (high efficiency due to the large available surface) and heterogeneous catalysis (easy recycling of the catalyst).

As is obvious, various applications of ferrofluids all require different properties and pose various constraints on the dispersions. To understand ferrofluid behaviour and to improve applications or develop new ones, careful study relating to fluid stability, control of surfactants, particle sizes, materials and physical behaviour are essential. In addition, ferrofluids, consisting of permanent spherical dipoles freely rotating in a nonmagnetic matrix can be regarded as a model for molecular dipoles, that are much more difficult to study.

9 References

- [1] Rosensweig R E 1985 *Ferrohydrodynamics* (Cambridge: Cambridge University Press).
- [2] Berkovski B, Bashtovoy, V. 1996 *Magnetic fluids and applications handbook* (New York: Begel House Inc.).
- [3] Charles S W and Popplewell J 1982 *Endeavour, New Series* **6** 153.
- [4] Lubbe A S, Alexiou, C., Bergemann, C. 2001 *J. Surgical Research* **95** 200.
- [5] Rosensweig R E 1982 *Scientific American* **247**(4), 136.
- [6] Geraedts, Wieringa (Océ) *personal communication*
- [7] Sun S, Murray C B, Weller D, Folks L and Moser A 2000 *Science* **287** 1989.
- [8] Stamm C, Marty F, Vaterlaus A, Weich V, Egger S, Maier U, Ramsperger U, Fuhrmann H and Pescia D 1998 *Science* **282** 449.

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Patricia Berger et al., „*Preparation and Properties of an Aqueous Ferrofluid*“, J. Chem. Educ. **76** (7), 943-948 (1999)

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A) Herstellung eines Fe_3O_4 - Ferrofluids (hydrophil)²

$\text{FeCl}_2 \cdot 4\text{H}_2\text{O}$ (3.29 g, 16.5 mmol) and $\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$ (8.68 g, 32.1 mmol) were dissolved together in 380 mL demineralised water. **Magnetite** (Fe_3O_4) was precipitated by adding 20 mL 25% NH_3 while stirring vigorously. After collection of the precipitate using a small permanent magnet, the supernatant was removed. The precipitate was washed with demineralised water. -> **Fe_3O_4 -FFaq**

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C) Herstellung eines $\gamma\text{-Fe}_2\text{O}_3$ - Ferrofluids (hydrophil)²

To the product of A) (**Fe_3O_4 -FFaq**) 40 mL 2 M HNO_3 (13,7 mL 65%ige HNO_3 auf 100 mL verdünnen) was added, and the mixture was stirred for 5 minutes. The oxidation to **maghemite** ($\gamma\text{-Fe}_2\text{O}_3$) was completed by adding 60 mL 0.35 M $\text{Fe}(\text{NO}_3)_3$ (14,14g $\text{Fe}(\text{NO}_3)_3 \cdot 9\text{H}_2\text{O}$ in 100 mL Lösung) to the mixture and stirring at its boiling point for one hour.

Next, the reddish sediment was washed with 2 M HNO_3 to keep the pH at a low value, after which it was redispersed in demineralised water. The resulting dispersion was coded **Fe_2O_3 -FFaq**.

D) Herstellung eines $\gamma\text{-Fe}_2\text{O}_3$ - Ferrofluids (lipophil)²

To be able to disperse particles in organic solvents, they were coated with oleic acid, using the method Ewijk²

2 mL of **Fe_2O_3 -FFaq** was flocculated by adding a few drops 25 % NH_3 . The sediment, collected by using a small magnet, was washed with demineralised water, after which approximately 10 mL water was added. Under mild stirring, 2 mL oleic acid was added, which caused the magnetic material to transfer to the oil. The black oil phase, separated from the water phase, was washed a few times with acetone and dried under nitrogen flow.

The resulting maghemite particles, coated with an oleic layer, were redispersed in cyclohexane. The resulting dispersion was coded **$\gamma\text{-Fe}_2\text{O}_3$ -FFch**.

² Ewijk van G A, Vroege G J and Philipse A P 1999 *J. Magn. Magn. Mater.* **201** 31-33